


For more than 20 years, Innovative Eyewear has been the pioneer in the field of color light therapy, and now we are proud to present our newest innovation:

### **PRiSMA SpektroChrom Color glasses!**

The SpektroChrom method for color light therapy has been developed by Dinshah P. Ghadiali more than 100 years ago, and is the best documented modern color light method.



This system, that attributes special colors to every organ and physiological function, has convinced us so much, that we optimized and fine-tuned our colors to suit the SpektroChrom method. From now on, we will only offer SpektroChrom color glasses.

We would like to invite you to get to know this method for yourself, and warmly recommend the website [www.spektrochrom.de](http://www.spektrochrom.de). Here you can get an insight into the SpektroChrom method, and you can get detailed treatment plans for more than 300 medical conditions, after a quick and free registration.

For all those who have used our “Vitalbrillen” before, and are now interested in SpektroChrom, there is good news: Since all our colors have always been chosen for optimum frequencies and efficacy, most of our colors are compatible with the SpektroChrom method. The missing colors can be bought in our [Upgrade Sets](#).

#### **Here are the changes in detail:**

SpektroChrom Yellow (Y) is a little more saturated than our old yellow.  
 The former Violett is called Magenta (M) in the SpektroChrom system.  
 SpektroChrom Purple (P) and Violett (V) are new colors.  
 The former Magenta is Scarlet (S) in the SpektroChrom system.

To upgrade a set of 11 Vitalbrillen to a full SpektroChrom set, you need 3 new colors:

SpektroChrom Yellow (Y), Violett (V), Purple (P).

An old 9 color set can be upgraded with 4 new colors:

SpektroChrom Yellow (Y), Violett (V), Purple (P), Lemon (L).